

Configuración básica FACS ARIA IIu

Láser 407 nm		
PMT	A	B
LP mirror	502LP	
BP filter	BP510/50	450/50
Fluorocromos	Alexa 430 CFP BD Horizon V500 Brilliant Violet 510 LIVE/DEAD Fixable Aqua	Cascade blue Pacific blue DAPI Hoeschst Alexa Fluor 405 eFluor 450 BD Horizon V450 SYTOX blue Cell Trace Violet stain

Láser 488 nm					
PMT	A	B	C	D	E
LP mirror	735LP	655LP	595LP	556LP	502LP
BP filter	BP780/60	BP695/40	BP610/20	BP575/26	BP530/30
Fluorocromos	PE-Cy7 PerCP eFluor 710	PerCP-Cy5.5 Ioduro de propidio	PE-Texas Red PE-CF594	PE Ioduro de propidio	FITC Alexa 488 GFP CFSE Fluo 3 Fluo 4 YFP SYTOX green LIVE/DEAD Fixable Green

Láser 633 nm		
PMT	A	B
LP mirror	735LP	
BP filter	BP780/60	BP660/20
Fluorocromos	APC-Cy7 APC-Alexa 750 Alexa 750	APC Alexa 633 Alexa 647 TO-PRO-3 Cy5 SYTOX Red FxCycle Far Red LIVE/DEAD Fixable Far Red